



PERFECTSTYLE

DERMALSTYLE®

MESOSTYLE®







Discover the Gentle Secret of Eternal Youth!

Nowadays, people are getting older and older. However, a youthful-appearing skin remains an ideal of beauty, and it contributes to a sense of well-being and self-confidence of many people. Like Theodor W. Adomo already said: "If you look good, you sell better."

With its different products – ranging from hyaluronic fillers, PDO threads for wrinkle correction, mesotherapy products to advanced skincare products – W&O medical esthetics GmbH, based in Oberursel, Germany, offers a safe, little invasive, and high-quality possibility to restore the long-term youthful and fresh appearance of the skin.

The root

The history of W&O medical esthetics GmbH starts with the special friendship of the two company founders: Dr. Najib Chichakli and Lais Mrad-Agua grow up together in Hama, Syria. Several years later, the friends' paths of life divide due to the politically difficult situation of the country. Dr. Najib Chichakli pursues medical studies in Austria; Lais Mrad-Agua pursues business studies in Germany. But though their lives take different courses, the friends always keep in touch. Years later, the perfect synergy of their professional qualifications and their unique mutual trust induce them to establish a company together.

The presen

For several years already, W&O medical esthetics GmbH has successfully distributed Dermalstyle® Fillers on international markets. With Perfectstyle (PDO-Threads) and Mesostyle® (Mesotherapy-Products) further products for wrinkle treatment and volume boost of the skin come in addition. And that's not the end: Dr. Najib Chichakli and Lais Mrad-Agua are still concentrating on revealing the secrets of eternal youth with new methods and products.

Renowned within the field of cosmetic treatment and with more than eight years of experience, Dr. Najib Chichakli and Lais Mrad-Agua and their team lead the way in anti-aging wrinkle correction and skin volume restoration procedures.

Taking the knowledge and experience Dr. Chichakli and his business partner and longtime friend Lais Mrad-Agua have now developed a line of advanced skincare products. The products are initially formulated as a cosmeceutical skin care range that complements advanced cosmetic dermatology treatments.

Although a steady growth of sales figures can be registered, Dr. Najib Chichakli and Lais Mrad-Agua attach great importance to the retention of their successful business strategy: a small business with excellent personal service as well as long-term business relations based on mutual trust.

On the following pages you will find all necessary information on our product range



Esthetical Fillers



Dermalstyle® Light contains the highest concentration of non-animal sourced and non-crosslinked hyaluronic acid currently available. Dermalstyle® Light is a viscoelastic solution to compensate the loss of hyaluronic acid due the ageing. It is especially designed for injection into the superficial dermal tissue, to improve tone and elasticity of the skin and/or to act as filler for small lines such as crow's feet, smile lines or smoke lines surrounding the mouth.



(

Concentration 18 mg/ml hyaluronic acid
Stabilizing agent 20 mg/ml glycerol
Dynamic viscosity [mPa.s] ~ 150.000
Molecular weight ~3 mio. Daltons
Other ingredients citrate buffer

Indications superficial wrinkles on face, neck, décolleté and back of hands Injection depth superficial to deep dermis

Technique multi-puncture (nappage), micro-papular technique **Needle** 30 G

Duration at the beginning 3 sessions (3 week interval), 2 to 3 month

Monaco Broschüre_RZ02.indd 2-3

Dermalstyle® Regular

Dermalstyle® Regular contains non-animal sourced, crosslinked hyaluronic acid and has a concentration of

to severe facial wrinkles and folds (such as nasolabial folds). It is recommended for all areas of the face Dermalstyle® Regular is indicated for injections into the mid to deep dermis for the correction of moderate where the wrinkles are not too pronounced.



Concentration 23 mg/ml hyaluronic acid

Reticulating agent BDDE

Dynamic viscosity [mPa.s] ~ 500.000

Molecular weight ~2,7 mio. Daltons

Other ingredients phosphate buffer, NaCl Indications moderately pronounced wrinkles in all facial areas

Technique linear, Criss-cross technique

Injection depth mid to deep dermis

Duration 6 - 9 month

Dermalstyle® Forte

Dermalstyle® Forte consists of non-animal sourced, crosslinked hyaluronic acid and with a concentration

subcutis or supraperiosteally. Dermalstyle® Forte is an injectable biodegradable gel to increase or restore the volume of the face, remodel facial contours and correct deeper wrinkles and folds. It is indicated for injections into the deep dermis,



(

Dynamic viscosity [mPa.s] ~1.000.000 Reticulating agent BDDE Other ingredients phosphate buffer, NaCl Molecular weight ~2,7 mio. Daltons Injection depth deep dermis Concentration 23 mg/ml hyaluronic acid **Indications** highly pronounced wrinkles Technique fanning, linear, depot technique

Duration 6 - 9 month





Dermalstyle® Forte Lidocaine

concentration of 23 mg/ml and 0,3% lidocaine. Dermalstyle® Forte Lidocaine consists of non-animal sourced, crosslinked hyaluronic acid and with a

deep dermis, subcutis or supraperiosteally. The additive anaesthetics enable a largely painless treatment. face, remodel facial contours and correct deeper wrinkles and folds. It is indicated for injections into the Dermalstyle® Forte Lidocaine is an injectable biodegradable gel to increase or restore the volume of the



Reticulating agent BDDE Concentration 23mg/ml hyaluronic acid and 0,3% lidocaine

Dynamic viscosity [mPa.s] ~1.000.000 Other ingredients phosphate buffer, NaCl Molecular weight ~2,7 mio. Daltons **Indications** highly pronounced wrinkles

Injection depth deep dermis Technique fanning, linear, depot technique

Duration 6 - 9 month Needle 27 G

Dermalstyle® Smile Lidocaine

a concentration of 23 mg/ml and 0,3% lidocaine. Dermalstyle® Smile Lidocaine contains a special non-animal sourced, crosslinked hyaluronic acid and has

correction of lip asymmetries. The additive lidocaine enables a largely painless treatment of the sensitive lip area and provides for more treatment comfort. Dermalstyle® Smile Lidocaine is especially indicated for volume augmentation, modelling, contouring and



(

Dynamic viscosity [mPa.s] ~550.000 Reticulating agent BDDE Molecular weight ~2,7 mio. Daltons Other ingredients phosphate buffer, NaCl Injection depth deep dermis Concentration 23 mg/ml hyaluronic acid and 0,3% lidocaine Indications lip contour and volume Technique serial puncture, linear, tenting technique Duration 6 - 9 month Needle 27 G



Monaco Broschüre_RZ02.indd 8-9

Advantages of Dermalstyle®

Dermalstyle® fillers are generally used to rejuvenate the skin, to treat facial wrinkles and folds, and to enhance the appearance of the lips.

Dermalstyle® products (Dermalstyle® Regular, Forte and Dermalstyle® Forte Lidocaine, Smile Lidocaine) belong to the safest dermal fillers on the market, due to an extremely low level of the crosslinking agent BDDE.

Dermalstyle® is a clear, non-animal-based, biodegradable gel that is composed of hypoallergenic hyaluronic acid. You do not need to test before treatment.

Dermalstyle® products are sterile, viscoelastic, clear, isotonic, biodegradable, and homogenized injectable gel implants.

They are highly versatile products that provide optimal aesthetic results and can be used for maximum correction in many different areas of the face such as nasolabial folds (smile lines), vertical lip lines, marionette lines (lip corners) and crow's feet.

Dermalstyle® offers a complete range of wrinkle-filling products. It consists of four different variations which can restore the facial volume, redefine and augment lips, repair cutaneous depressions, and revitalize the dermis.

(

The innovative unique production process gives Dermalstyle® exceptional viscosity and elasticity; thus it integrates homogeneously into the skin for a smooth, natural, revitalized appearance.

The benefits of Dermalstyle®:

- no allergy testing required
- stimulates collagen synthesis
- maintains skin hydration
- increases skin's suppleness & firmness
- long-lasting effect
- natural youthful look

Reflecting the latest state of research, the ingredients of Dermalstyle® are extremely pure, biocompatible and non-allergenic. Thus our products significantly over-comply with the stringent European Directive for medical devices (CE-certification) standards. Clinical studies confirm the optimum compatibility, safety and efficiency of Dermalstyle®:


Quality with long-lasting effects.

PERFECTSTYLE

PDO Threads

Double impact for an effective result

Perfectstyle uses 100% absorbable threads made of polydi-oxanone (PDO) which has already been in use in surgery for decades. The threads are differently structured – smooth, spiral, or with microscopically small hooks.



(

These special dissolving threads are applied to the hypoder-mic fat tissue. Thus the skin is tightened and the tissue is lifted to its original position.

The treatment with biodegradable polydioxanone threads combines two effects: the immediately visible tightening of the treatment area is followed by a reactivation – so-called bio-stimulation – of the connective tissue so that the skin wins additional firmness. Systematic wrinkle treatment with Perfectstyle can thus lead to an immediate, naturally appearing and long-lasting skin lifting.

Monaco Broschüre_RZ02.indd 10-11

MONO GRIP

- Straight thread with a sharp 35° needle
- Universally suitable for face, neck, cleavage, and body
- Smoothing of small wrinkles; enhancement of skin volume, skin tightening
 Suitable for all skin types, especially for thinner skin



SINGLE SCREW

- Straight, spirally twisted thread with a sharp 35° needle
- Universally suitable for face, neck, cleavage, and body
- Smoothing of medium-sized wrinkles, stronger enhancement of skin volume and skin tightening compared to Mono Grip threads
- Suitable for normal skin.



DOUBLE SCREW

(

- 2 Screw threads twisted to one thread, with a sharp 35° needle
- Universally suitable for face, neck, cleavage, and body
- Smoothing of distinct, deep wrinkles; strongest enhancement of skin volume and most intense skin tightening compared to Mono Grip and Single Screw
- Suitable for normal and thicker skin



COG MASTER GEAR

- PDO thread with bi-directional barbs in the form of small cones ("spear heads") with new I-blunt needle
- Universally suitable for face, neck, and body
- Suitable for all skin types, especially for very thin skin with little subcutaneous
- fatty tissue and in case of sagging fascia
- Collagen synthesis and immediately visible lifting effect

COG MASTER GEAR





COG Pro1/COG Pro2/COG Pro3

- Multi-directional 360° COG threads in 3-D design
- Suitable for skin, neck, and body
- Smoothing of wrinkles and enhancement of skin volume as well as strong lifting and pulling effect

COG PRO I (Sharp Cannula)



















Monaco Broschüre_RZ02.indd 12-13

COGGING SCREW with barbs

- $\bullet~$ Spiral PDO threads with bi-directional barbs and a sharp, rounded 45° needle
- Suitable for face, neck, and body
- Suitable for all skin types
- Collagen synthesis and immediately visible lifting effect

COGGING SCREW



COGGING SCREW PLUS

Needle inside Cog with gripping and retention direction vector






Visu


(

of the needle interior view

- Short, strong PDO threads with barbs and new I-blunt needle
- Suitable for body and face, especially for nose corrections
- Also suitable for the enhancement of skin volume and thus for smoothing deeper wrinkles

Visu I







for the treatment of cellulite, hyperpigmentation, hair loss and scars.

With Mesostyle® products you can guarantee your customers immediately visible, natural and effective

To achieve the desired outcome, the right choice and combination of high-quality substances is essential

(

skin look younger, fresher, tighter and more vitalised. The targeted insertion of nutrients across the skin bar-

rier allows the process to take full effect in the desired areas. Therefore, this therapy is especially effective

Esthetical mesotherapy actively supports cell regeneration as well as the production of collagen, and it

stops the moisture loss of the skin by rehydrating it with substances such as hyaluronic acid. This makes the

Even without a scalpel, a significant rejuvenation and optimisation of the tissue structures can be achieved.

Modern esthetical dermatology offers various ways to improve the complexion and the look of the skin.

Mesostyle[®] – *Safe, Gentle And Effective*

Esthetical mesotherapy – the boost for your skin

MESOSTYLE®

Mesostyle® products are perfectly suited for skin rejuvenation of face, neck, décolleté and hands.

MESOSTYLE® Radiance C

Skin regeneration and tightening (face/body)

AREAS OF APPLICATION

- Inhomogeneous skin colouration
- Problems with hyperpigmentation
- Dry and tired skin
- Light-damaged skin

INGREDIENCES:

Gluconate. Aqua, Ascorbic Acid, Sodium Hyaluronate, Manganese Gluconate, Glycerin, Superoxide Dismutase, Sodium



MESOSTYLE® Lipo X

Effective reduction of smaller fat depots and cellulite (face/body)

AREAS OF APPLICATION:

- Fibrous cellulite
- Localised fat depots

INGREDIENCES:

Officinalis, Rutin, Taurin, Sodium Gluconate. Aqua, Carnitin, Theophylline, Cynara Scolymus, Melilotus



MESOSTYLE® Hair Plus

Treatment of thin, weak hair and hair loss

AREAS OF APPLICATION:

- Dehydrated hair fibres
- Slow hair growth

INGREDIENCES:

Aminophylline, Cumarin, Sodium Gluconate. phosphate, Pyridoxine Hcl, Panthenol, Biotin, Taurine, Aqua, Zinc Gluconate, Thiamine Hcl, Adenosine Tri-



MESOSTYLE® Stretch Control

For more skin elasticity and firmness (face/body)

AREAS OF APPLICATION:

• Decreased firmness, atony of the skin of face and body

Loose skin and flaccid cellulite

Stretch marks

INGREDIENCES:

Niacinamide, Sodium Gluconate. Acid, Methylsilanol Mannuronate, Copper Gluconate, Aqua, Sodium Hyaluronate , Meclofenoxate, Thioctic



MESOSTYLE® HA PURE 2%

For more hydration and revitalization (face and body)

AREAS OF APPLICATION:

- Deyhydrated skin
- Photoaged skin

INGREDIENCES:

Aqua, Sodium Hyaluronate, Sodium Gluconate.

Wrinkled and sagging skin



Who can be treated with MESOSTYLE® products?

(

ingredient complex is individually composed, thus providing the right product for each indication. The active substance complexes of MESOSTYLE® products can be used for a variety of treatments as each

What are the benefits of a MESOSTYLE® therapy?

- No inactive periods due to treatment
- Quick, lasting and sustained effects
- Effective results
- No preparation time
- How many treatments are necessary?

of treshness. on, in 14-day intervals. It is, however, also possible to receive a single treatment now and then for a boost We recommend a course of at least three treatments, depending on the requirements and the skin indicati-

Monaco Broschüre_RZ02.indd 14-15

MESOSTYLE® Microneedling-System

Auto-Microneedling-System ULTIMA A6

treatment of scars. Stimulation of the collagen synthesis for the smoothing of wrinkles, the rejuvenation of the skin and the

The so-called micro-needling method is versatile and suitable for the following indications:


- smoothing and tightening of the skin
- skin rejuvenation, smoothing of wrinkles, treatment of sun-damaged skin
- any kind of scars (especially due to burns and acne)
- stretch marks
- hyper- and hypopigmentation
- reduction of pore sizes and refinement of the complexion

treatment, new sterile, sharp needles are used, which considerably reduces pain during treatment. it with a depth between 0.2 and 2.5 mm, according to the respective skin depth. For every new course of tion. It contains of disposable caps fitted with microneedles that pierce the skin vertically and penetrate The ULTIMA A6 microneedling system stimulates the body's own physiological abilities of skin regenera-

Benefits of the auto micro-needling system Ultima A 6

- Prevention of infections due to disposable needle caps
- Needle length adjustable between 0.2 and 2.5 mm
- Wireless, rechargeable model with 2 lithium batteries
- 5 speed levels: level 1: up to 3900 revolutions per minute; level 5: up to 15000 revolutions per minute

• To be used with both batteries and mains supply



DERMALSTYLE® COSMETICS

Your daily skin care

Daily regimenfor skin ageing prevention

technology, Dermalstyle® skincare offers a powerful regime for the skin. synergy. Formulated from a multitude of specially selected ingredients using innovative advanced skincare The sophisticated technology of our formulas include multiple ingredients which work in balance and Every Dermalstyle® product represents the most advanced and most effective skincare technology available.

HYDROLIFT Eye Serum

Intensely hydrating and perfecting eye serum

puffiness. Hydrolift Eye Serum provides a radiantly fresh and natural look expression lines. It helps to diminish the appearance of dark circles under the eyes and nic acid. Peptides provide regeneration and help to smooth the signs of wrinkles and effective ingredient complex helps to stimulate the skin's natural synthesis of hyaluro-Hydrolift Eye Serum is an innovative moisturizing serum with triple effect. The highly

> STYLE' Hydrolift Eye Serum

> > **(**

- Extremely hydrating
- Regnerating and revitalizing
- Helps to reduce small wrinkles



Intense care for a visibly relaxed and regenerated skin complexion.

the skin against harmful stress and environmental factors natural regeneration processes and its moisture balance. Furthermore it helps to protect signs of premature aging. Additional anti-aging ingredients help to optimize the skin's time it helps to prevent the formation of new wrinkles and protects the skin against soothe the skin and to diminish the appearance of wrinkles and fine lines. At the same The innovative and highly effective serum, with an exclusive peptide complex, helps to

- Helps to soften the appeareance of wrinkles and lines
- Promotes the moisture barrier
- Helps to improve skin structure



Monaco Broschüre_RZ02.indd 16-17

-

HYDROLIFT Creme

Restores a youthful natural glow to tired and dry skin

helps to protect the skin against harmful environmental influences. facial contours and supports skin regeneration. In addition, a repair complex visible lifting effect. A special ingredient complex ensures an improvement of This skin-optimizing moisturizer for demanding skin was developed for an

- Visible Lifting Effect
- Extremely hydrating
- Regenerating effect



HYDRODERM Repair Creme

(

Soothing care and effective support for skin in need of regeneration

skin's natural moisture barrier is strengthened. effect of post anti-wrinkle injections is supported by this topical cream and the relaxed skin feeling and supports the natural regeneration process. The visible especially after small medical aesthetics treatments. It provides a positive and This post-treatment cream is specially adapted to the needs of stressed skin,

- Soothing advantages
- Regenerates the skin
- Extremely hydrating



collagen synthesis and at the same time increases skin elasticity. redness. The combination of natural active agents, consisting of hyalu quickly soothes and relieves the skin in case of irritations, tautness an cells and reduces moisture loss. The innovative formula of the mask natural hygroscopic agent hyaluronic acid stimulates and revives the for all skin types soothes the skin and intensively moisturises it. The ronic acid, aloe vera, minerals and vitamins (A, C and E), activates the

The Dermalstyle® Cooling, Calming & Rehydration Mask

- Immediate moisture balancing with cooling effect
- Reduces swelling
- Protects against free radicals
- Smoothens and strengthens the skin



The Dermalstyle® Regeneration & Revitalisation Mask

(

of peptides, vitamins (A, C and E) and herbal stem cells provide the skin against UV light and free radicals. skin with highly effective anti-ageing agents. Antioxidants protect the type. High and low molecular weight hyaluronic acids, a combination is a moisturizing, vitalizing and regenerating face mask for every skin

- Regenerating, revitalising
- Improves the elasticity and tonicity of the skin
- Protects against free radicals
- Reduces small wrinkles



Monaco Broschüre_RZ02.indd 18-19



08.03.18 19:27